



SARS-CoV-2 Rapid Antigen Test

Testing ready to go when you need to know

SARS-CoV-2 Rapid Antigen Test allows for decentralised testing at the point of care and assists in expanding the range and volume of virus testing into diverse settings, particularly if laboratory testing is limited or unavailable. In addition, the SARS-CoV-2 Rapid Antigen Test is an instrument-free test, enabling testing in many settings including rural and low-infrastructure areas.

The SARS-CoV-2 Rapid Antigen Test enables fast decision making to determine if patients need to be put in quarantine, thus reducing the risk of further spreading of the virus. It also allows for screening of individuals after confirmed exposure to a SARS-CoV-2 infected person or individuals that are high risk for exposure such as health care workers.

Main benefits

- » **No analyser** required
- » **Fast results** with a 15 minute turn around time
- » **Prefilled reagents** to ensure simple ease of use, minimal handling steps, and reduced risk of manual errors
- » **Virus inactivation** is achieved in the buffer after 2 minutes
- » **Sample options:** test sample directly from the patient or to test sample already contained in viral or universal transport media
- » **Storage:** test kits can be at room temperature

Getting fast answers is just a few steps away

Performing the Roche SARS-CoV-2 Rapid Antigen Test

1

Collect the nasopharyngeal sample with the supplied swab.

2

Insert the swab into an extraction buffer tube and mix.

3

Remove the swab while squeezing the sides of the tube and replace cap.

4

Apply 3 drops of extracted sample to the test device.

4

15 – 30 min

Risk of incorrect results. Do not read the test result after 30 min.

Read the test result at 15 - 30 min.

Delivering reliable results at the point of care

Clinical Performance

Clinical performance was determined from a sample cohort of 975 samples across two prospective studies. The overall relative sensitivity was 95.5% (Ct value ≤ 30) and specificity 99.2%,^{1,2}

SARS-CoV-2 Rapid Antigen Test

Sensitivity	95.5% (Ct value ≤ 30 ; 95% CI: 91.8% - 97.8%)
Specificity	99.2% (95% CI: 98.2% - 99.7%)


Reference:

1. Compared to an FDA EUA authorised laboratory based real time reverse transcription PCR test.
2. SARS-CoV-2 Rapid Antigen Test Product Information AU V1.0

Assay Characteristics

Test type	Qualitative - rapid chromatographic immunoassay
Sample type	Nasopharyngeal <ul style="list-style-type: none"> ▪ Direct collection from patient ▪ Sample contained in VTM/UTM
Target antigen	Nucleocapsid (N)
Time to result	15 minutes
Time to read result	From 15 - 30 minutes
Storage temperature	2 - 30°C
Stability (test, opened pouch)	1 hour once the test has been opened
Registration	Included on ARTG 352250

Understanding the results when you need to know



Test Description	Quantity Per Kit	Catalogue Number
SARS-CoV-2 Rapid Antigen Test	25	09327592190
SARS-CoV-2 Rapid Antigen Control Kit	10 sets of positive and negative	09338322190



All trademarks mentioned herein are protected by law.
© 2021 Roche Diagnostics

Roche Diagnostics Australia Pty Limited
ABN 29 003 001 205
2 Julius Avenue
North Ryde NSW 2113 Australia
Phone: 02 9860 2222
E-mail: australia.dia_cardiology@roche.com.au
MC-AU-00937

COVID-19 rapid antigen point of care tests should be conducted by a medical practitioner or a registered/enrolled nurse who can provide an individual with appropriate advice and treatment.